



Bure Standard

The Bure Standard is one of the markets original combined sit-to-stand and rehabilitation aids, with a design refined over two decades. Eminently suitable within both the acute and post acute environments, the Bure Standard creates improved ergonomic working conditions for care givers.

The Bure Standard is available in electric, hydraulic and manual height adjustment, with all models providing stepless adjustment of the arm support.

BURE STANDARD

Ergonomic handles

Ergonomically designed handles that are steplessly adjustable in every direction for optimum individual comfort. This is extremely important for stroke patients who can thus reach their "usual" hand position.

Height-adjustable handles

See accessories.

Sustainable

StandTall Walkers are designed to last. For example, they have extra reinforcement bands around the arm rest cushions to prevent them from splitting from impacts to the side.

Stable

The walker has anatomical pads that provide relief for shoulders, arms and the neck and shift the focus to larger muscle groups such as leg and chest muscles instead. This improves stability and makes the walker easier to steer.

Easily reachable hand controller

The hand controller is located such that both of the carer's or user's hands can easily reach it.

Flexible stepless adjustment

Electrically driven height adjustment makes it easy to pre-set the StandTall Walker to the desired height. Once the patient is standing in position the height can be steplessly adjusted for comfort and the proper support for an optimal walking configuration.

Power Rise

Bure Standard Electric provides excellent assistance when helping patients and users stand. The strain on personnel, and thus repetitive strain injuries, are minimised. Some patients can even get up on their own with the assistance of the StandTall Walker.

Brakes and directional locks

There are optional one or two-handed brakes that can be positioned so that the user can reach them from the walking position. The handbrake acts on the rear castors and the brake assembly also includes directional locks that make it easier for the user to stay on course when walking.

Open design

The open design provides the same stability but with greater room for manoeuvre. The patient does not feel shut in and it makes assistance with standing up and dressing easier.

Waterproofing

All Bure Standard Electric parts have IP 65 ratings and can therefore be used in wet spaces.



Electric

Bure Standard electric is ideal for users who need a walker equipped with power rise and which are capable of height adjustment once a patient is standing.



Hydraulic

Bure Standard gas is ideal for users who need a quality walker capable of stepless height adjustment but who do not need assistance from a power rise actuator.



Manual

Bure Standard manual is great for users who need good, individually adaptable walkers but who do not need to adjust the height very often and who have no need of assistance from a power rise actuator.

Technical Specification

Part No.	Product	Width	Max user weight	Length	Height	Armrest/Internal dimension
56-312L-US	Bure Standard, Electric, 125mm Castors	78cm	150kg	80cm	90-130cm	25-50cm
56-311	Bure Standard, Hydraulic, 125mm Castors	78cm	150kg	80cm	90-129cm	25-50cm
56-310-KD	Bure Standard, Manual, 125mm Castors	78cm	150kg	80cm	90-133cm	25-50cm

All walkers may be ordered with 75, 100 or 125mm castors – specify when ordering. (Unless otherwise specified, the heights shown below refer to products with 125 mm castors).