



TR 1000/-W

Hygiene Chair

*Hygiene Equipment
Designed for People*

TR 1000/-W

The TR 1000/-W battery operated hygiene chair is height adjustable and can be tilted backward and forward for easy and efficient transfers and showering.



Registered at Medical Products Agency, Sweden as a MDR Class 1 product according to MDR 2017/745.

TR reserves the right to make design changes and/or improvements without notice.



The TR 1000/-W height adjustable hygiene chair is a great solution to minimise static load and offers a safe ergonomic working environment for caregivers.

The full power feature of the TR 1000/-W allows the caregiver to use the hand control to adjust the height and tilting of the chair, so that they can be in a good position when giving care to the individual. This provides a safe patient handling environment for both the resident and the caregiver.





TR 1000/-W

Hygiene Chair

Hygiene Equipment
Designed for People

- Clean design of chassis and chair makes access easy during transfers and showering. It also makes the TR 1000/-W easy to clean.
- Suitable to use with most sling hoists and standing/raising aids to make transfers as smooth and safe as possible.
- Seat and chassis create possibilities for toileting over standard toilet or with a bedpan attached to the chair.
- The rear chair design features a smooth surface for ease of cleaning.
- Belt over hip feature for patient safety.
- The adjustable and foldable footrest ensures security and comfort for the residents (available option).



Bedpan, comes with a lid.



Height adjustable and foldable footrest.

Specifications

Lifting stroke length	400 mm / 16"
Working load	200 kg / 440 lbs
Tilt backwards	25 degrees
Tilt forward	7 degrees
Structure	Powdercoated, zinc plated steel and stainless steel.
Castors	Twin 100 mm / 4" locking mechanism
Control	Hand control, for height adjustment and tilt.
Battery	Removable and rechargeable LI-Ion 24V/4Ah

Models:

0521000	TR1000 Shower chair, straight armrests
0521010	TR1000 Shower chair, front closing armrests
0511000	TR1000-W Shower chair, straight armrests extended width (50mm)
0511010	TR1000-W Shower chair, front closing armrests (50mm)

Options/Accessories:

0521015	Adjustable and foldable footrest
0521020	Chest belt
850169	Bed pan including lid
3232622	Wallbracket for battery pack

