



TR 2810

Shower panel

Hygiene Equipment
Designed for People

TR 2810 Complete shower system

TR 2810 is a shower panel with thermostatic mixing valve and a separate cleaning system that together with TR shower trolleys fulfill a complete shower system.

This system means care is comfortable, safe and hygienic for both patients and caregivers.

The cleaning system runs on cold water to minimise the risk of inhalation when cleaning equipment.

The shower panel is also available with a mechanical lock for cleaning detergent (option).



Registered at Medical Products Agency, Sweden as a MDR Class 1 product according to MDR 2017/745.

TR reserves the right to make design changes and/or improvements without notice.





TR 2810

Shower panel

Hygiene Equipment
Designed for People

- Easily operated thanks to clear and clean design
- Integrated hygienic and safe cleaning system for easy cleaning of equipment
- Built in bottle for cleaning detergent which is automatically mixed with cold water in the cleaning system
- Thermometer on the shower hose is available as option to read the water temperature.
- Hot and cold water supply for patient shower with thermostatic mixer, valve for patient shower, shower spray handle and white shower hose.
- Cold water supply for cleaning with safety valve for equipment cleaning, cleaning spray handle, yellow cleaning hose shower, integrated cleaning, cleaning flow valve setting inside panel.

Specifications

Width	580 mm
Height	360 mm
Depth	130 mm
Recommended mounting from floor	1250 mm



Display water temperature

