



TR 2810

Shower Panel

*Hygiene Equipment
Designed for People*

Spare part list



Registered at Medical Products
Agency, Sweden,
as a MD Class 1 product,
MDD 93/42EEC.

TR reserves the right to make
design changes and/or
improvements without notice.



Versions produced from June 2009

Please specify serial no. and article no. when contacting TR Equipment

TR 2810 Art. no. 0520070-02




TR EQUIPMENT



TR 2810

Shower Panel

Hygiene Equipment
Designed for People



Head office

TR Equipment AB, Sweden
Phone: +46 140 38 50 60
info@trequipment.se
www.trequipment.com

Group Subsidiaries

TR Group Inc, US, Phone: +1 (800) 752 6900
TR Equipment Ltd, UK, Phone: +44 (0) 844 335 8386
TRanspoRehaMed GmbH, Germany, Phone: +49 (0) 212 645850-0



TR EQUIPMENT



TR 2810

Shower Panel

*Hygiene Equipment
Designed for People*

When ordering spare parts please specify model and batch number, see label on the shower panel.

Position	Art No	Description
1	829284	Thermostatic mixer
2	829288	Valve shower
3	829289	Cleaning valve incl. key
4	829307	Valve stem 1/2"
5	829289-02	Key for cleaning valve
10	829171	Shower handle
11	829172	Shower handle bracket
12	829097	Shower hose white
13	829096	Cleaning hose yellow
14	850614	Cleaning bottle
15	7426293	Lid, cleaning bottle
16	829069	Watts non return valve
17	7230166	Cover plastic
	829169	ORAS Handshower on/off white
	829284-02	Cartridge assembly R215
18	829070	Cleaning mixture setting injektor

Head office

TR Equipment AB, Sweden
Phone: +46 140 38 50 60
info@trequipment.se
www.trequipment.com

Group Subsidiaries

TR Group Inc, US, Phone: +1 (800) 752 6900
TR Equipment Ltd, UK, Phone: +44 (0) 844 335 8386
TRanspoRehaMed GmbH, Germany, Phone: +49 (0) 212 645850-0



TR EQUIPMENT