



TR Stretcher for TR Comfortline

Hygiene Equipment
Designed for People

Dressing and shower stretcher for TR Comfortline

The shower stretcher is especially designed for placing above a Comfortline bathtub to be used as showerbench or changing table.

For increased safety and security, the shower stretchers can be equipped with safety rails on one or both sides. The shower stretcher is easy to put on and take off and is made from strong and durable materials that are easy to keep clean.



Registered at Medical Products Agency, Sweden as a MDR Class 1 product according to MDR 2017/745.

TR reserves the right to make design changes and/or improvements without notice.



Specifications

Length: 1 665 mm

Width: 795 mm

Weight: 8,5 kg

Weight, siderail: 3 kg per siderail

Working Load Limit: 150 kg

Design: Flame retardent PVC and acrylic coated polyester fabric, on a stainless steel powder coated frame



TR EQUIPMENT



TR Stretcher for TR Comfortline

Hygiene Equipment
Designed for People

Features

- The shower stretcher is easy to put on and take off from the bathtub
- It is made from strong and durable materials that are easy to keep clean
- The fabric is stretched tightly and firmly around the frame
- The stretcher is design with smooth shapes and surfaces providing for easy maintenance and good hygiene

Options

- Safety Side Rails, for increased safety and security. The side rails can be easily mounted on one or both sides of the stretcher. Safety rails can be folded all the way down into a vertical position.



Art. no.	Model
05184000	TR dressing and shower stretcher
05184001	TR Stretcher, right siderail
05184002	TR Stretcher, left siderail
05184003	TR Stretcher, complete with 2 siderails

