



TR Swing bathtub

Reclining and height adjustable

Hygiene Equipment
Designed for People

TR Swing Reclining bath with side access

The Swing is a free standing, reclining bathtub designed for a satisfying and relaxed bathing experience.

The TR Swing can be equipped with several Wellness options, such as LED-lighting, Air Spa- and bluetooth music system.

To make it easier for the care giver, the TR Swing can be equipped with Height adjustment, Autofill function and Cleaning system.

Simple installation of the TR Swing allows it to be positioned in almost any central bathing areas. Controls have been designed with simplicity in mind and the smooth design makes it easy to clean and maintain.



Registered at Medical Products Agency, Sweden as a MDR Class 1 product according to MDR 2017/745.

TR reserves the right to make design changes and/or improvements without notice.





TR Swing bathtub

Reclining and height adjustable

Hygiene Equipment
Designed for People

- Innovative and attractive design
- Easy to operate and learn, quick introduction
- Compact design for easy placement
- Ergonomic tub shape
- Easy to clean and maintain
- Detached, no wall or floor fixing
- Hand shower with long hose
- Overflow control and pop-up drainage including draining hose
- White as standard color
- Two grab handles inside the tub
- Reinforced water hoses for water installation

Specifications

| | |
|------------------------|-----------------------------------|
| Length | 2040 mm |
| Width | 800 mm |
| Access height | 570 mm |
| Lifting speed | 10 mm/s / 0,3"/s |
| Working load limit | 200 kg |
| Power supply | 230 V, 50 Hz, 10 A / 110 V, 60 Hz |
| Drainage outlet | 40 mm |
| Overall weight (empty) | 140 kg |

Models:

TR Swing reclining bath

TR Swing reclining and height adjustable bath

Options/Accessories:

- Autofill
- Air spa with 12 nozzles
- Bluetooth/Aux music system with integrated speakers
- Integrated Cleaning system
- Led light with button control



Reclined bath



Bluetooth and AUX music system

